



Shenzhen Zhongchuang 3D Technology Co., Ltd

## Technical Data Sheet (TDS)--OPPOSITE JAW MODEL Resin

Version: V1.0

Issue Date: 2023/3/15

Product: OPPOSITE JAW MODEL Resin

### 1. Company Information

|                        |                                                                                                                                 |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Company Name</b>    | Shenzhen Zhongchuang 3D Technology Co., Ltd.                                                                                    |
| <b>Company Address</b> | Room 204, Building 1, Jiatai Industrial Park, No. 1 Industrial Zone, Tangwei Community, Fuhai Street, Bao'an District, Shenzhen |
| <b>Brand Name</b>      | ZonMed                                                                                                                          |
| <b>Email</b>           | info@zonmedtech.com                                                                                                             |
| <b>Website</b>         | www.zonmedtech.com                                                                                                              |

### 2. Product Information

|                          |                                                                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name</b>      | Repair-Model Resin                                                                                                                                             |
| <b>Product Model</b>     | PRS-R-19                                                                                                                                                       |
| <b>Application Scope</b> | Used for printing upper and lower opposing jaw models, suitable for basic diagnostic scenarios such as occlusal relationship analysis and denture arrangement. |

### 3. Composition on Ingredients

| Substance name                                                      | Conc. (%)   | CAS.       |
|---------------------------------------------------------------------|-------------|------------|
| (Octahydro-4,7-methano-1H-indenediyl)bis(methylene) bismethacrylate | 20.00-45.00 | 43048-08-4 |
| Diurethane dimethacrylate, mixture of isomers                       | 20.00-35.00 | 72869-86-4 |
| Poly(ethyleneglycol)dimethacrylate)                                 | 10.00-30.00 | 25852-47-5 |
| Polyoxyethylene diacrylate                                          | 10.00-30.00 | 26570-48-9 |
| 1-Hydroxycyclohexyl phenyl ketone                                   | 2.00-5.00   | 947-19-3   |

### 4. Physical & Chemical Properties

| Test Item       | Standard Value               | Test Method           |
|-----------------|------------------------------|-----------------------|
| Appearance      | Beige                        | Visual Inspection     |
| Viscosity (25℃) | 2200-5300CPS                 | Rotational Viscometer |
| Density (25℃)   | 1.13 ~ 1.19g/cm <sup>3</sup> | Density Meter         |
| pH Value        | 6.5 ~ 7.5                    | pH Meter              |
| Curing Mode     | Light-curing                 | Operational Test      |
| Curing Time     | 2.5s                         | Curing Machine        |

## 5. Performance Parameters (After Curing)

| Parameter                   | Standard   | Unit    | Result    |
|-----------------------------|------------|---------|-----------|
| Shore                       | ASTM D2240 | Shore D | 72        |
| Bending strength            | ASTM D790  | Mpa     | 75-85     |
| Bending modulus             | ASTM D790  | Mpa     | 1950-2115 |
| Notch impact strength       | ASTM D256  | Mpa     | 38-40.5   |
| Breaking elongation rate    | ASTM D638  | %       | 7-10.2    |
| Water absorption rate       | ASTM D570  | %       | 0.8-0.98  |
| Hot deformation temperature | /          | ° C     | <73       |

## 6. Usage Instructions

- ◆ Preparation: Clean and dry the dental treatment area before use to ensure no moisture and oil pollution.
- ◆ Print Pre-processing: Thoroughly shake the resin to eliminate air bubbles and pour it into a dental light-curing printer (385 – 405 nm). Add supports to the model, configure the appropriate print parameters, level the build platform, and start printing.
- ◆ Part Removal and Cleaning: Remove the printed part and ultrasonically clean it in 95% isopropyl alcohol for 2 – 3 minutes; dry thoroughly with a cool air stream.
- ◆ Curing: Place the part in a UV curing machine, ensuring exposure from multiple angles for a total duration of at least 3 – 5 minutes; standard: not feeling sticky when holding in hand.
- ◆ Grinding and Polishing: Remove supports, adjust occlusion and margins, and polish in stages to achieve a natural tooth-like luster.
- ◆ Seal the bottle tightly after each use to prevent liquid volatilization and contamination.

## 7. Storage & Transportation Conditions

- ◆ **Storage Temperature:** 5°C - 25°C, sealed storage
- ◆ **Shelf Life:** 24 months (unopened and stored in standard conditions)
- ◆ **Transportation Requirements:** Avoid extrusion, collision, high temperature and sunlight exposure
- ◆ **Storage Notice:** Keep away from fire, heat sources and strong oxidants

## 8. Package Specification

OPPOSITE JAW MODEL: 1kg/bottle



## **9. Safety Description**

- ◆ This product is dedicated to dental clinical use, forbidden for oral ingestion and external use on non-dental parts.
- ◆ Avoid contact with eyes. In case of accidental contact, rinse with plenty of clean water immediately.
- ◆ If allergic reaction occurs during use, stop using immediately and take corresponding medical measures.
- ◆ Waste materials after use shall be disposed of in accordance with medical waste management regulations.

## **10. Disclaimer**

All data in this Technical Data Sheet is based on our current test results and industry standards, which is for reference only. The purchaser shall judge the applicability of the product independently according to the actual use environment. Our company shall not be liable for any accidental loss caused by improper operation and non-compliant use. We reserve the right to update and modify the TDS data without prior notice.