

DLP Light-curing 3D Printer

FEATURES

- 1 Printing accuracy $\pm 40.5\mu\text{m}$
- 2 Imported industrial DMD optical machine
- 3 Pixel level micro-adjustment and anti-aliasing function
- 4 Printing speed up to 30 minutes
- 5 Z-axis accuracy up to 0.0035mm



DEVICE PARAMETERS

Molding size	144*81*172mm
Resolution	1920*1080
Layer thickness	0.02-0.2mm
Molding speed	10-30mm/h
Screen size	7 inches
Light source	DMD chip; Imported 405nm light source
Equipment power	40W
Software functions	Timing lock and offline unlocking, Automatic support function, Modify thickness angle function

DEVICE APPLICATION



Full Mouth Movable Base Resin



Temp C&B Resin




Permanent Crown Resin



Flexible Base Resin



Clear Aligner Resin



Surgical Guide Resin