



Ceramix Ceramic Stain Full Kit

Product Datasheet | Brand: ZonMed | Origin: Shenzhen, China | Certification: CE, ISO 13485

1. Basic Information

Product Name	Ceramix Ceramic Stain Full Kit (G01-G20 Series 20 Bottles)
Model	PRS-S-G
Category	Dental Stain and Glaze
Series	Ceramix Stain
Contents	18 standard colors G01-G18 + G19 Glaze (FP) + G20 Glaze (20 bottles)

2. Selling Point Keywords

- Ceramix G01-G20 Full Kit
- Complete Ceramic Stain Set
- Zirconia Colorant Set
- Glass Ceramic & Press Ceramic Colorant Set

3. SEO Elements (English, for overseas marketing)

SEO Title	Ceramix Ceramic Stain Full Kit G01 G20 Low Fusing Paste for Dental Zirconia
SEO Keywords	Ceramix Stain Kit, G01 G20 Ceramic, Zirconia Colorant
SEO Description	Ceramix Ceramic Stain full kit G01 to G20 colors zirconia and glass ceramic. Single bottles also available. Contact us for wholesale pricing.

4. Basic Attributes

Brand	ZonMed
Model	PRS-S-G
Origin	Shenzhen, China
Certification	CE, ISO 13485

5. Ordering & Trade Attributes

MOQ	1 Kit
Price	Negotiable
Payment	T/T in advance (payment before shipment)
Supply Ability	10,000 bottles / month
Delivery	3-7 working days
Packaging	Carton, 230 mm x 270 mm x 70 mm



6. Custom Attributes

Product Type	Ceramic Stain Kit
Series	Ceramix Stain
Color Code Range	G01-G20
Number of Colors	20
Net Content	5 ml / bottle
Sintering Temperature	Zirconia 800C; Glass Ceramic / Press Ceramic 750-760C
Heating Rate	50C / minute
Total Cycle	Approx. 15 minutes

7. Product Description

Ceramix Ceramic Stain Full Kit - Complete G01-G20 Series

The Ceramix Ceramic Stain Full Kit contains the full series of 20 standard colors with unlimited free mixing, enabling dental laboratories to reproduce any natural zirconia / glass ceramic / press ceramic shade with a stable high-gloss finish.

Kit Contents (20 bottles)

- 18 standard colors G01-G18: Bright Red, Dark Red, Pink, Light Pink, Sepia, Bronze, Khaki, Orange, Yellow, Green, Dark Green, Dark Blue, Light Blue, Lilac, Black, Gray, Beige, White (fluorescent)
- G19 Glaze (FP) (fluorescent)
- G20 Glaze (non-fluorescent)

Core Advantages

- Universal: compatible with all ceramic materials and staining techniques
- Flexible: various combinations mix to create any shade
- Easy to use: innovative paste formula for hassle-free application
- Realistic: stable high-gloss finish with natural mimicry

Sintering Guide

- Zirconia: 800C; Glass ceramic (press and CAD): 750-760C
- Heating rate: 50C per minute; total cycle approx. 15 minutes

Single Bottle Purchase

Single bottles (individual shades) are available for purchase. To order any individual color from G01 to G20, please contact us.

Application

Zirconia crowns and bridges, glass ceramic (press / CAD) restorations, veneers and other ceramic restorations.

Packed in a color box, 230 mm x 270 mm x 70 mm. Wholesale pricing available - contact us.

8. Product Images

