



Shenzhen Zhongchuang 3D Technology Co., Ltd.

## PolyFix Polymer Stain Full Kit

Product Datasheet | Brand: ZonMed | Origin: Shenzhen, China | Certification: CE, ISO 13485

### 1. Basic Information

|              |                                                                                                |
|--------------|------------------------------------------------------------------------------------------------|
| Product Name | PolyFix Polymer Stain Full Kit (R01-R21 Series   21 Bottles)                                   |
| Model        | PRS-S-R                                                                                        |
| Category     | PolyFix Stain                                                                                  |
| Series       | PolyFix Stain                                                                                  |
| Contents     | 18 standard colors R01-R18 + R19 Diluent + R20 Transparent Glaze + R21 Soft Glaze (21 bottles) |

### 2. Selling Point Keywords

- PolyFix R01-R21 Full Kit
- Complete Resin Stain Set
- Denture Base Colorant Set
- PMMA Colorant Set
- PEEK Colorant Set
- Polymer Colorant Set

### 3. SEO Elements (English, for overseas marketing)

|                 |                                                                                                                                          |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------|
| SEO Title       | PolyFix Polymer Stain Full Kit R01 R21 Light Cure Resin Colorant for Denture Base                                                        |
| SEO Keywords    | PolyFix Stain Kit, R01 R21 Resin Stain, Denture Base Colorant                                                                            |
| SEO Description | PolyFix Polymer Stain full kit R01 to R21 for denture base resin. Single bottles also available. Contact us for wholesale pricing today. |

### 4. Basic Attributes

|               |                 |
|---------------|-----------------|
| Brand         | ZonMed          |
| Model         | PRS-S-R         |
| Origin        | Shenzhen, China |
| Certification | CE, ISO 13485   |

### 5. Ordering & Trade Attributes

|                |                                          |
|----------------|------------------------------------------|
| MOQ            | 1 Kit                                    |
| Price          | Negotiable                               |
| Payment        | T/T in advance (payment before shipment) |
| Supply Ability | 10,000 bottles / month                   |



|                  |                                    |
|------------------|------------------------------------|
| <b>Delivery</b>  | 3-7 working days                   |
| <b>Packaging</b> | Color box, 230 mm x 270 mm x 70 mm |

## 6. Custom Attributes

|                         |                                                    |
|-------------------------|----------------------------------------------------|
| <b>Product Type</b>     | Light-Cure Resin Stain Kit                         |
| <b>Series</b>           | PolyFix Stain                                      |
| <b>Color Code Range</b> | R01-R21                                            |
| <b>Number of Colors</b> | 21                                                 |
| <b>Net Content</b>      | 3 ml / bottle                                      |
| <b>Curing Time</b>      | 30 seconds or less                                 |
| <b>Application</b>      | Denture base, temporary crown & bridge, PMMA, PEEK |
| <b>Shelf Life</b>       | 24 months                                          |

## 7. Product Description

### PolyFix Polymer Stain Full Kit - Complete R01-R21 Series

The PolyFix Polymer Stain Full Kit achieves natural appearance restoration and personalized aesthetic customization of removable denture bases, temporary bridges and crowns through precise color matching and creative color shaping.

### Kit Contents (21 bottles)

- 18 standard colors R01-R18: Bright Red, Dark Red, Pink, Light Pink, Sepia, Bronze, Khaki, Orange, Yellow, Green, Dark Green, Dark Blue, Light Blue, Lilac, Black, Gray, Beige, White
- R19 Diluent
- R20 Transparent Glaze
- R21 Soft Glaze

### Key Features

- High color stability with intelligent color matching technology
- Fast curing speed of 30 seconds or less
- Multi-shade combination to create any custom shade
- Antibacterial, wear-resistant and corrosion-resistant

### Single Bottle Purchase

Single bottles (individual shades) are available for purchase. To order any individual color from R01 to R21, please contact us.

### Application

Gum, full mouth movable base, PMMA block, PEEK and other polymers.

Packed in a color box, 230 mm x 270 mm x 70 mm. Wholesale pricing and OEM support available - contact us.

[illegible]