



Sintering Paste

Product Datasheet | Brand: ZonMed | Origin: Shenzhen, China | PRS Ecosystem

1. Basic Information

Product Name	Sintering Paste
Model	PRS-P-S01
Series	PRS Ecosystem - Auxiliary Material
Net Content	14 g

2. Selling Point Keywords

- Sintering Paste Fixation
- Dental Zirconia Support Paste
- Ceramic Crown Firing Paste

3. SEO Elements (English, for overseas marketing)

SEO Title	Dental Sintering Paste Prevent Discoloration for Zirconia Ceramic Crown Firing Process
SEO Keywords	Sintering Paste, Zirconia Firing Paste, Dental Crown Support
SEO Description	Dental Sintering Paste fixes zirconia and ceramic crowns during firing with no residue. Prevents discoloration and saves drying time. Contact us now.

4. Basic Attributes

Brand	ZonMed
Model	PRS-P-S01
Origin	Shenzhen, China
Certification	CE / ISO 13485 (to be confirmed)

5. Ordering & Trade Attributes

MOQ	1 box
Price	Negotiable
Payment	T/T in advance (payment before shipment)
Supply Ability	1,000 boxes / month
Delivery	3-7 working days
Packaging	Standard packaging



6. Product Description

Sintering Paste is suitable for veneers, inlays, onlays, partial crowns, single crowns, bridges and other restorations that are inconvenient to place directly on the sintering tray or sintering pin during porcelain baking or glazing. It helps stabilize the restoration, keep the shade stable and save operation time.

Product Features

- **Stable Fixing:** secures hard-to-place restorations during porcelain furnace sintering
- **Prevents Discoloration:** restorations stay color-fast after firing, no yellowing or graying
- **Easy Removal:** peels off easily after sintering, no sandblasting, zero residue
- **Saves Drying Time:** no extra pre-drying needed in the porcelain furnace

Application

Porcelain baking / glazing of veneers, inlays, onlays, partial crowns, single crowns, bridges and other ceramic restorations.

Net content: 14 g. Note: The above is factory information; please refer to the actual product in case of upgrades.

7. Product Images

