



DENTAL 3D PRINTING RESIN

Provide a complete set of specialized
dental solutions you need

TEMPORARY C&B RESIN **REPAIR MODEL RESIN** **NIGHT GUARD PAD RESIN** **CLEAR ALIGNER RESIN** **WATER WASHABLE MODEL RESIN**
TRAY RESIN **SURGICAL GUIDE RESIN** **OPPOSITE JAW MODEL RESIN** **FLEXIBLE BASE RESIN** **IBT RESIN** **WAX CASTING RESIN**
TRY-IN RESIN **FULL MOUTH MOVABLE BASE RESIN** **BIOMIMETIC GUM RESIN** **PLANTING MODEL RESIN**
TEMPORARY VENEER RESIN **SMILE VENEER RESIN** **ORTHODONTIC MODEL RESIN** **PERMANENT CROWN RESIN**

<https://www.zonmedtech.com>
<https://www.chinadental3dprintingresin.com>

Zhongchuang Medical Group Co., Ltd

compatible with



ASIGA



rapidshape



ELEGOO



ANYCUBIC



RAYSHAPE




UNIFORMATION

and the other open system printers...

Founded in 2019, ZonMed is a leading high-tech enterprise dedicated to the research, development, and production of the full PRS Ecosystem.

ZonMed PRS™ provides a seamless, closed-loop digital workflow from initial scan to final restoration.

With our own advanced ensure strict quality our innovation. We market directly, reliable and perior solutions to



manufacturing facility, we control and continu-serve the global exporting our aesthetically su-dental profession-

and regions. ZonMed your trustable partner for restorative materials that deliver consistent results worldwide.

PRS™ ECOSYSTEM

is committed to being high-performance

5000+

Cooperative Labs & Clinics

2019

Establishment

95%

Customer Retention Rate



DIGITAL INFRASTRUCTURE

Precise 3D printing systems and smart post-processing units optimized for stable, high-speed dental production.



HIGH-PRECISION RESINS

Medical-grade resins featuring an ultra-precise marginal fit of < 50µm for crowns and surgical guides.



AESTHETIC FINISHING

Rapid low-temperature glaze and staining systems designed to replicate natural tooth optics in minutes.



ESSENTIALS

High-quality finishing tools, polishing kits, and clinical accessories to complete the PRS ecosystem.


Catalogue




PRS-R-01 PERMANENT CROWN 01



PRS-R-02 TEMPORARY C&B 02



PRS-R-03 TRY-IN 03



PRS-R-04 FLEXIBLE BASE 04




PRS-R-05 FULL MOUTH MOVABLE BASE 05




PRS-R-06 TEMPORARY VENEER 06



PRS-R-07 SMILE VENEER 07



PRS-R-08 WAX CASTING 08



PRS-R-09 IBT 09





PRS-R-10 TRAY

10



PRS-R-11 SURGICAL GUIDE

11



PRS-R-12 BIOMIMETIC GUM

12




PRS-R-13 NIGHT GUARD PAD

13




PRS-R-14 CLEAR ALIGNER

14



PRS-R-15 WATER WASHABLE MODEL

15



PRS-R-16 OPPOSITE JAW MODEL

16



PRS-R-17 ORTHODONTIC MODEL

17



PRS-R-18 PLANTING MODEL

18



PERMANENT CROWN

No.: PRS-R-01



PRODUCT FEATURES

Porcelain Powder Reinforced Formula

High-Precision Detail Restoration

Multi-Layer Aesthetic Effect

Strong Fracture Resistance




Used for fabricating single crowns and inlays, with high strength, wear resistance, and natural aesthetics, replacing traditional metal-porcelain restorations.

Specification 1:
0.25kg/bottle(contains 60% ceramic powder)

Color: BL, A1, A2, A3, A3.5, B1, C1

Specification 2:
0.5kg,1kg/bottle (contains 10% ceramic powder)



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	90-96
Bending strength	ASTM D790	Mpa	180-208
Bending modulus(Specification 2)	ASTM D790	Mpa	4300-4700
Bending modulus(Specification 1)	ASTM D790	Mpa	7000-8800
Density(Specification 1)	/	g/cm ³	1.6-1.7
Notch impact strength	ASTM D256	Mpa	18-28
Breaking elongation rate	ASTM D638	%	4-6.5
Water absorption rate	ASTM D570	%	0.35-0.42
Tensile strength	ASTM D256	Mpa	85-90

TEMPORARY C&B

No.: PRS-R-02

PRODUCT FEATURES

- High-Strength Temporary Support
- Natural Light-Transmitting Aesthetics
- Strong Stability



Used for temporary crown/bridge restoration after root canal treatment or tooth preparation to protect abutment teeth and maintain chewing function until permanent restorations are completed.



Color: BL, A1, A2, A3, A3.5, B1, C1

Specification: 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	86-90
Bending strength	ASTM D790	Mpa	148
Bending modulus	ASTM D790	Mpa	3200-3600
Notch impact strength	ASTM D256	Mpa	20-40
Breaking elongation rate	ASTM D638	%	8-11
Water absorption rate	ASTM D570	%	0.6-0.95
Tensile strength	ASTM D790	Mpa	80-95

TRY-IN

No.: PRS-R-03



PRODUCT FEATURES

- Efficient Trial Fitting
- Convenient On-Site Adjustment



Producing disposable trial dental bridges, allowing patients to preview restoration outcomes and assisting dentists in optimizing occlusion.

- Color: A2
- Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	78
Bending strength	ASTM D790	Mpa	101-115
Bending modulus	ASTM D790	Mpa	2450-2650
Notch impact strength	ASTM D256	Mpa	30.5-38.5
Breaking elongation rate	ASTM D638	%	5.5-6.5
Water absorption rate	ASTM D570	%	0.12-0.19
Hot deformation temperature	/	℃	<65

FLEXIBLE BASE

No.: PRS-R-04

PRODUCT FEATURES

Superior Elasticity

Gum-Friendly

Low Hardness

Super-Elastic Clasp



Specifically designed for patients with sensitive gums or compromised alveolar bone, our flexible base evenly distributes forces to minimize discomfort.



Color: Pink

Specification: 0.5kg, 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	69-73
Bending strength	ASTM D790	Mpa	35-45
Bending modulus	ASTM D790	Mpa	600-790
Tensile strength	ASTM D256	Mpa	30-35
Breaking elongation rate	ASTM D638	%	75-95
Water absorption rate	ASTM D570	%	0.4-0.52

FULL MOUTH MOVABLE BASE

No.: PRS-R-05



PRODUCT FEATURES

High Strength

Long-Lasting Performance

Precise Force Distribution



Premium complete denture bases for a perfect fit-delivering stability, restored function, and a natural look for confident smiles.

Color: Pink

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	86
Bending strength	ASTM D790	Mpa	142
Bending modulus	ASTM D790	Mpa	2650-2850
Notch impact strength	ASTM D256	Mpa	40.5-58.5
Breaking elongation rate	ASTM D638	%	7.5-9.0
Water absorption rate	ASTM D570	%	0.6-0.82
Hot deformation temperature	/	°C	<70

TEMPORARY VENEER

No.: PRS-R-06

PRODUCT FEATURES

Ultra-Thin Design

High-Transparency Bright White

Efficient Molding



Used to make short-term trial or transitional veneers for patients to experience the restorative effects, assist dentists in adjusting aesthetic plans, or provide temporary protection for tooth whitening.



Color: BL

Specification: 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	92
Bending strength	ASTM D790	Mpa	168-188
Bending modulus	ASTM D790	Mpa	4800
Notch impact strength	ASTM D256	Mpa	39-55
Breaking elongation rate	ASTM D638	%	7.5-8.5
Water absorption rate	ASTM D570	%	0.5-0.75
Hot deformation temperature	/	℃	<70

SMILE VENEER

No.: PRS-R-07



PRODUCT FEATURES

Adjustable Comfortable Fit

Long-Lasting Stain Resistance

Natural Light-Transmitting Texture

Temperature-Sensitive Deformation, Portable



Directly print removable whitening veneers to quickly improve tooth color and shape without tooth grinding, suitable for daily beauty.

Color: BL

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	58
Bending strength	ASTM D790	Mpa	35-42
Bending modulus	ASTM D790	Mpa	850-950
Tensile strength	ASTM D256	Mpa	25.5
Breaking elongation rate	ASTM D638	%	92
Water absorption rate	ASTM D570	%	0.45-0.58
Hot deformation temperature	/	°C	<60

WAX CASTING


No.: PRS-R-08

PRODUCT FEATURES

- No Residue, Easy Casting
- Extreme Detail Restoration
- Rapid Demolding, Time-Saving



Used for printing wax patterns in metal framework/porcelain-fused-to-metal crown casting and ceramic casting, replacing traditional manual wax carving to improve efficiency.



- Color: Red
- Specification: 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Sintering temperature	/	℃	800-900
Bending strength	ASTM D790	Mpa	70-78
Storage time of finished products	/	H	<6
annealing temperature	/	/	1-1.5
Residual rate	/	%	<0.15

IBT

No.: PRS-R-09



PRODUCT FEATURES

Segmented Printing For High Adaptability

Precision Molding

Transparent, Invisible, And Easy To Clean



Used for fabricating fixed orthodontic brackets and accessories to assist in invisible aligner/traditional bracket orthodontics and accurately transmit orthodontic forces.

Color: Transparent

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	75
Bending strength	ASTM D790	Mpa	65-78
Bending modulus	ASTM D790	Mpa	1760-1950
Notch impact strength	ASTM D256	Mpa	23.5-28.5
Breaking elongation rate	ASTM D638	%	6.5-9.5
Water absorption rate	ASTM D570	%	0.65-0.89
Hot deformation temperature	/	°C	<75

TRAY


No.: PRS-R-10

PRODUCT FEATURES

- High Strength And Toughness
- Strong Material Compatibility
- Precise Positioning And Efficient Operation



Used for fabricating customized trays such as orthodontic trays and impression trays to meet the needs of precise impression taking and positioning in oral treatments.



- Color: Green
- Specification: 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	86
Bending strength	ASTM D790	Mpa	98-105
Bending modulus	ASTM D790	Mpa	2250-2750
Notch impact strength	ASTM D256	Mpa	28.5-35.5
Breaking elongation rate	ASTM D638	%	5.5-6.5
Water absorption rate	ASTM D570	%	0.13-0.19
Hot deformation temperature	/	°C	<60

SURGICAL GUIDE

No.: PRS-R-11



PRODUCT FEATURES

High Construction Precision

High Toughness And Anti-Fracture

Transparent, Visible, And High Precision



Used for printing dental implant/orthodontic surgical guides to precisely position surgical pathways and assist dentists in performing complex surgical procedures.

Color: Transparent

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	82
Bending strength	ASTM D790	Mpa	75-85
Bending modulus	ASTM D790	Mpa	1880-2050
Notch impact strength	ASTM D256	Mpa	25.5-35.5
Breaking elongation rate	ASTM D638	%	6.5-8.5
Water absorption rate	ASTM D570	%	0.5-0.85
Hot deformation temperature	/	°C	<75

BIOMIMETIC GUM

No.: PRS-R-12

PRODUCT FEATURES


High Elasticity and Tear Resistance

Excellent Dynamic Fit Performance

Natural Gingiva Aesthetics



Used for fabricating simulated gingiva, applied in dental implant surgery simulation and denture aesthetic design to restore the shape and elasticity of natural gingiva.



Color: Pink

Specification: 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	35-38
Bending strength	ASTM D790	Mpa	2.5-5.0
Bending modulus	ASTM D790	Mpa	1250-1750
Notch impact strength	ASTM D256	Mpa	0.8-1.0
Breaking elongation rate	ASTM D638	%	200-320
Water absorption rate	ASTM D570	%	0.13-0.19

NIGHT GUARD PAD

No.: PRS-R-13



Transparent Glaze (soft)

PRODUCT FEATURES

Heat-Activated Shaping

Effortless Wearing

Anti-Fracture Toughness

Corrosion-Resistant And Easy To Maintain



Used to fabricate nighttime bruxism guards and occlusal splints to cushion the impact of grinding, protecting teeth.

Color: Transparent blue

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	70-75
Bending strength	ASTM D790	Mpa	35.5-48
Bending modulus	ASTM D790	Mpa	1280-1490
Notch impact strength	ASTM D256	Mpa	45.5-48.5
Breaking elongation rate	ASTM D638	%	60-80
Water absorption rate	ASTM D570	%	0.40-0.65

CLEAR ALIGNER

No.: PRS-R-14

PRODUCT FEATURES

- Portable Printing Process
- Aesthetics With Effective Correction
- Precise Mechanical Force Transmission



Transparent Glaze (soft)

Directly print invisible aligners to achieve tooth movement through gradual force application, meeting the aesthetic needs of adult orthodontics.



- Color: Transparent
- Specification: 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	72-75
Bending strength	ASTM D790	Mpa	35-45
Bending modulus	ASTM D790	Mpa	600-800
Notch impact strength	ASTM D256	Mpa	25-35
Breaking elongation rate	ASTM D638	%	80-100
Water absorption rate	ASTM D570	%	0.5-0.7

WATER-WASHABLE MODEL

No.: PRS-R-15



PRODUCT FEATURES


- Direct Tap Water Cleaning
- Eco-Friendly Process
- Low-Cost & High Efficiency
- Multi-Scenario Application



For printing basic dental models, suitable for entry-level scenarios such as simple models.

Color: Gray

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	85.7
Bending strength	ASTM D638	Mpa	53
Bending modulus	ASTM D638	Mpa	1498
Tensile modulus	ASTM D638	Mpa	1355
Breaking elongation rate	ASTM D638	%	11.8
Water absorption rate	ASTM D570	%	0.82-0.95

OPPOSITE JAW MODEL

No.: PRS-R-16

PRODUCT FEATURES

Boil-Resistant & Crack-Free

Stable After Boiling


No Release Agent Required



For printing opposing dental models, used in basic diagnostic scenarios such as occlusal relationship analysis and denture arrangement.

Color: Beige

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	72
Bending strength	ASTM D790	Mpa	75-85
Bending modulus	ASTM D790	Mpa	1950-2115
Notch impact strength	ASTM D256	Mpa	38-40.5
Breaking elongation rate	ASTM D638	%	7-10.2
Water absorption rate	ASTM D570	%	0.8-0.98
Hot deformation temperature	/	°C	<73

ORTHODONTIC MODEL

No.: PRS-R-17



PRODUCT FEATURES

High-Temperature Resistance
For Thermoforming

Strong Film Compatibility

Ultra-Low Viscosity & High Fluidity



For printing orthodontic correction models, supporting orthodontic processes such as in visible aligner thermoforming and bracket positioning.

Color: White

Specification: 1kg/bottle



PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	85
Bending strength	ASTM D790	Mpa	80-89
Bending modulus	ASTM D790	Mpa	1650-1850
Notch impact strength	ASTM D256	Mpa	15-20.5
Breaking elongation rate	ASTM D638	%	9-11.2
Water absorption rate	ASTM D570	%	1.0-1.18
Hot deformation temperature	/	°C	<180

PLANTING MODEL


No.: PRS-R-18

PRODUCT FEATURES

- High-Precision Detail Replication
- Excellent Stability
- Strong Impact Resistance
- Compatible With Sectioned Die Models



For printing high-precision dental implant surgical models, used for modeling complex anatomical structures.



- Color: Yellow, Gray
- Specification: 1kg/bottle

PRODUCT PERFORMANCE

Parameter	Standard	Unit	Result
Shore	ASTM D2240	Shore D	80
Bending strength	ASTM D790	Mpa	80-90
Bending modulus	ASTM D790	Mpa	2050-2315
Notch impact strength	ASTM D256	Mpa	36-41.5
Breaking elongation rate	ASTM D638	%	7-10
Water absorption rate	ASTM D570	%	0.68-0.92
Hot deformation temperature	/	℃	<80

INNOVATIVE RESINS

PERFECT SMILES



chongxiangli155@gmail.com

<https://www.zonmedtech.com>

<https://www.chinadental3dprintingresin.com>

B10,RM 2703 NEWTREND CENTRE,704PRINCEEDWARD ROAD
EAST,SAN POKONG, KLN,HK

Room 204, Administration Building, Jia Tai Pharmaceutical, No.1 Industrial Zone,
Tangwei Community, Fuhai Street, Bao'an District, Shenzhen City



Wechat



WhatsApp

